

Running the Updated HIREXSR Analysis Code (9/16/10)

The following should be added to your .bashrc file

```
LD_LIBRARY_PATH=$LD_LIBRARY_PATH:/home/labombard/idl_lib/fortran
```

To compile the code run @hirexsr_ini.bat and compile *widgets* as necessary (use .compile *widgetname*).

- 1) run **HIREXSR_UPDATE_TREE** which will create the HIREXSR node in spectroscopy and the ANALYSIS, CALIB and INFO structures.

- 2) fill in INFO data from .info files using **HIREXSR_LOAD_INFO2TREE** or copy from another shot using **HIREXSR_COPY_INFO**

MOD #1 _____ MOD #2 _____ MOD #3 _____ MOD #4 _____

- 3) write rest wavelength table to new shot using **HIREXSR_WRITE_WAVELENGTHS**

- 4) a) PERFORM CALIBRATIONS USING LOCKED MODE

1. compile and run the calibration widget *w_hirexsr_calib* for each module.

MOD #1 _____ MOD #2 _____ MOD #3 _____ MOD #4 _____

2. verify or modify the crystal detector alignment using the alignment widget

w_hirexsr_det_align

MOD #1 _____ MOD #2 _____ MOD #3 _____ MOD #4 _____

3. write white field calibration **HIREXSR_WRITE_WHITE**

MOD #1 _____ MOD #2 _____ MOD #3 _____ MOD #4 _____

4. write filter transmission **HIREXSR_WRITE_TRANS**

- b) COPY CALIBRATION DATA from another shot using **HIREXSR_COPY_CALIB**

- 5) write module ordering using **HIREXSR_WRITE_MORDER**

- 6) write binning using **HIREXSR_WRITE_BINNING** or copy from another shot using **HIREXSR_COPY_BINNING**

HE-LIKE _____ H-LIKE _____

- 7) run **HIREXSR_CALC_DATA** to form averages spectra and moments, write them to the tree and use *w_hirexsr_he_moments* or *w_hirexsr_h_moments* to examine the spectra and moments

W(0) _____ X(1) _____ Z(2) _____ LYA1(3) _____

- 8) run **HIREXSR_INVERT_DATA** to calculate profiles, write them to the tree and use *w_hirexsr_profiles* to examine the inversions.

W(0) _____ X(1) _____ Z(2) _____ LYA1(3) _____